

OPEN BASE, POWER, CHANGING TABLE



Only Clinton has the manufacturing know-how and the focus to engineer and bring to market a changing table that meets needs and saves money. Building on our reputation for dependable and affordable medical equipment, Clinton's Open Base, Power, Changing Table is easy to operate and comfortable for patients. It features sturdy steel construction, smooth electric power and two padded side rails that easily lock into place for safety.

Why choose anything else, when practical performance comes at such an affordable price?

- Easy-clean, one-piece top
- Two lock-in-place, padded side rails
- Four, smooth-roll, locking, 3" casters
- Electric power height adjustment with foot control
- Ultra-strong, uni-frame base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Seamless, rounded corner, upholstered top
- 2" (5 cm) ultra-firm padding
- Four, locking casters
- Two Lengths: 72" and 60"
- Made in the USA by Clinton Industries
- 400 lbs. (181 kg) load capacity under normal use



087-6 Extra-Height Option

- Height extension to accommodate Hoyer Lift
- Adds 6.25" to the height of the table from top pad to floor



80110-60

Overall

Length: 60"

Width: 27"

Height Range: 18.5" – 35.5"



80110-72

Overall

Length: 72"

Width: 27"

Height Range: 18.5" – 35.5"



Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black



CLINTON
INDUSTRIES INC.

800-441-9131 • 717-848-3519

Toll-Free Fax 866-522-1233 • Fax 717-843-5871

www.clinton-ind.com

California residents please be advised of the following, pursuant to Proposition 65:

⚠ WARNING: Products featured on this page can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

CLI-MA-079